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## REFERENTIES

**Adcock DM, Kressin DC, Marlar RA.**

Effect of 3.2% vs 3.8% sodium citrate concentration on routine coagulation testing. *Am J Clin Pathol* 1997; 107: 105-110.

**Adcock DM, Kressin DC, Marlar RA.**

The effect of time and temperature variables on routine coagulation tests. *Blood Coagulation and Fibrinolysis* 1998; 9: 463-470.

**Adcock DM, Kressin DC, Marlar RA.**

Minimum specimen volume requirements for routine coagulation testing. Dependence on citrate concentration. *Am J Clin Pathol* 1998; 109:595-599.

**Brien WF, Schaus MR, Cooper KE, O'Keefe BT, Inwoon M.**

Lupus anticoagulant testing: effect of platelet count on the activated partial thromboplastin time. *Br J Biomed Sc* 1993; 50: 114-116.

**Brandt JT, Triplett DA, Alving B, Scharer I.**

Criteria for the diagnosis of lupus anticoagulants: an update. On behalf of the Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibody of the Scientific and Standardisation Committee of the ISTH. *Thromb Haemost* 1995; 74:1185-1190.

**Capel P, Chatelain B, Leclercq R, Lust A, Masure R, Arnout J.**

Quality control in haemostasis. *Acta Clinica Belgica* 1992; 47: 308-318.

**Contant G, Gouault-Heilmann M, Martinoli JL.**

Heparin inactivation during blood storage: its prevention by blood collection in citric acid, theophylline, adenosine, dipyridamole –C.T.A.D. mixture-. *Thromb Res* 1983; 31:365-374.

**European Committee for clinical laboratory standards. Standard for specimen collection.** Part 2: Blood specimen by venipuncture. *ECCLS Document vol. 4, no. 1 (1987)*. Beuth Verlag GmbH, Berlin.

**Fiebig EW, Ezzell JE, Ng VL.**

Clinically relevant differences in prothrombin time and INR values related to blood sample collection in plastic vs glass tubes. *Am J Clin Pathol*. 2005;124:902-9.

**Gilmer PR.**

Preanalytical variables in coagulation testing. In: DA Triplett (editor), Standardization of coagulation assays: an overview. *College of American Pathologists. Skokie, Illinois, 1982; pp 1-8*.

**Gogstad GO, Dale S, Brosstad F, Brandsnes Ø, Holtlund J, Mørk E, Gärtner E, Borch SM.**

Assay of D-dimer based on immunofiltration and staining with gold colloids. *Clin Chem* 1993 ; 39 :2070-2076.

**Greaves M, Cohen H, Machin SJ, Mackie I.**

Guidelines on the investigation and management of the antiphospholipid syndrome. *Br J Haematol* 2000; 109: 704-715.

**Ho CH, Wu S-Y.**

The influence of time, temperature and packed cell on activated partial thromboplastin time and prothrombin time. *Thromb Res* 1991; 62: 625-633.

**Hope E, Mayorga SR, Robinson S, Goldberg L, Leveson JE, Marengo-Rowe AJ, Peerschke EIB.**

Preparation of plasma for coagulation testing. Evaluation of the StatSpin high-speed centrifuge. *Laboratory Medicine* 1991; 22: 190-193.

**Hillyard CJ, Blake AS, Wilson K, Rylatt DB, Miles S, Bunch R, Elms MJ, Barnes A, Bundesen PG.**

A latex agglutination assay for D dimer: evaluation and application to the diagnosis of thrombotic disease. *Clin Chem* 1987; 33:1837-1840.

**Iversen LH.**

Pre-analytical variation in the measurements of sensitive markers of coagulation and fibrinolysis: the influence of venipuncture and mixing of blood.

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*Infusionstherapie* 1984; 11: 46-50.

**Lippi G, Salvagno GL, Montagnana M, Guidi GC.**

Short-term venous stasis influences routine coagulation testing. *Blood Coagul Fibrinolysis* 2005; 16: 453-458.

**Lippi G, Salvagno GL, Poli G, Guidi GC.**

Influence of centrifugation temperature on routine coagulation testing. *Clin. Chem.* 2006; 52: 537-538.

**Lutze G, Breyer J, Naumann Ch, Zawta B**

Useful Facts about Coagulation, second edition, Roche Diagnostics GmbH, Mannheim 2000

**NCCLS**

**Procedures for the collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition** Number H3-A5. Pennsylvania, NCCLS; 2003; vol 23, No 32

**CLSI**

Collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemaostasis assays; Approved guideline – Fifth Edition. *CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute;; 2008.*

**Peterson P, Gottfried EL.**

The effects of inaccurate blood sample volume on prothrombin time (PT) and activated partial thromboplastin time (aPTT). *Thromb Haemost* 1982; 47:101-103.

**Reneke J, Ezzell J, Leslie S, Ng VL, Gottfried EL.**

Prolonged prothrombin time and activated partial thromboplastin time due to underfilled specimen tubes with 109 mmol/l (3.2%) citrate anticoagulant.

*Am J Clin Pathol* 1998; 109: 754-757.

**Ray MJ, Carroll PA, Just SJE, Hanson GAT.**

A low volume specimen container suitable for monitoring the aPTT of heparinized patients. *Blood Coagul Fibrinol* 1993; 4: 805-807.

**Sletnes KE, Gravem K, Wisloff F.**

Preparation of plasma for the detection of lupus anticoagulants and antiphospholipid antibodies. *Thrombosis Res* 1992; 66: 43-53.

**Stibbe J, Hemker HC, van Creveld S.**

The inactivation of factor VIII in vitro. *Thromb Diathes haemorrh* 1972; 27:43-58.

**Triplett DA.**

Screening for the lupus anticoagulant. *Res Clin Lab* 1989; 19: 379-389.

**Van den Besselaar AMHP, Chantarangkul V, Tripodi A.**

A comparison of two sodium citrate concentrations in two evacuated blood collection systems for prothrombin time and ISI determination. *Thromb Haemost* 2000 ; 84 :664-667.

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Effect of magnesium contamination in evacuated blood collection tubes on the prothrombin time test and ISI calibration using recombinant human thromboplastin and different types of coagulometer. *Thromb Res* 2005 ; 115 :239-244.

**Van Geest-Daalderop JHH, Mulder AB, Boonman-de Winter LJM, Hoekstra MMCL, Van den Besselaar AMHP.**

Preanalytical variables and off-site blood collection : influences on the results of the prothrombin time/international normalized ratio test and implications for monitoring of oral anticoagulant therapy. *Clin Chem* 2005 ; 51 :561-568.

**Walker ID.**

Blood collection and sample preparation: pre-analytical variation. In: J Jespersen, RM Bertina, F Haverkate (editors), 2<sup>nd</sup> revised edition of ECAT assay procedures. *A manual of laboratory techniques. Kluwer Academic Publishers, Dordrecht, 1999; pp 21-28*

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**Woodhams B, Girardot O, Blanco MJ, Colesse G, Gourmelin.**

Stability of coagulation proteins in frozen plasma. *Blood Coagul Fibrinolysis* 2001; 12:229-236.

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
$$0,20 - 0,13 = 0,07 \text{ ml citraat per ml plasma.}$$

Reagents

Select the curve for the total volume of anticoagulated blood required (10, 5, 2, or 1 mL). Enter the chart at the patient's packed cell volume on the horizontal axis and read off the corresponding volume of anticoagulant solution on the vertical axis. Place the volume in a collection tube and add blood up to the required total volume. The chart assumes a normal packed cell volume of 43%; in that case, 1 volume of anticoagulant is made up to 10 volumes with blood.

Alternatively, the anticoagulant volume may be calculated from the expression:

$$X = (100 - PCV) / (595 - PCV) \text{ vol,}$$

where, X is the volume of anticoagulant required to prepare unit volume of anticoagulated blood, and PCV is the packed cell volume in %; e.g., to determine the volume required for 5 mL anticoagulated blood, calculate 5 x mL.

(Reproduced from Ingram, G.I.C., Brozovic, M., and Slater, N.G.P. Bleeding disorders, investigations and management, 2nd edition. Blackwell's Scientific Publications, Oxford: pp. 244-245, 1982, by permission of the authors and publishers.)

CLSI H21-A5